

# Virosart® HC

## The High Capacity Virus Filter for Blood & Plasma Applications

### Product Information

Virosart® HC is especially developed for plasma derivatives in cooperation with one of the market leaders within this industry. The unique surface modified hydrophilic PES membrane can be easily wetted and shows a constant flow even with hydrophobic feed streams e.g. IVIG. Flexibility is given by using either cartridges in existing stainless steel housings or disposable capsules.

### Description

Choose your perfect fit from the Sartorius virus clearance strategy summarizing orthogonal technologies, manufacturing solutions, validation support and consultancy. The orthogonal technologies from Sartorius consist of virus inactivation as well as virus removal by chromatography and virus filtration.

The Virosart® product range includes four different virus retentive membranes, in order to provide the best solution for every application. Virosart® HC targets the removal of both small non-enveloped viruses (20 nm) e.g. PPV, MVM and larger enveloped viruses (> 50 nm) e.g. MuLV from a biopharmaceutical product, in particular from hydrophobic feed streams.



## Application & Positioning

The main applications for Virosart® HC for virus filtration are hydrophobic feed streams such as IVIG or any other plasma derivative as well as hydrophobic antibodies and recombinant proteins (< 150 kDa). Virosart® HC is used at the end of the purification process for virus filtration of the biopharmaceutical product.

At this stage the purity of the biopharmaceutical product is the highest and probability of virus filter blockage due to contaminants (DNA, CHOP, aggregates and lipoproteins) is the lowest. Although contaminants should be removed during the polishing process of the target molecule, small amounts might be sufficient to cause premature blockage of the final virus filter.

To prevent this, an efficient pre-filtration step, such as the Virosart® Max\*, might be required as protection for the Virosart® HC membrane. The optimum pre-filter to final filter ratio has to be identified during development of the process step.

## Product Benefits

Virosart® HC provides efficient virus removal to the biopharmaceutical product. Based on the unique surface modified double layer 20 nm PES membrane, Virosart® HC provides excellent capacity even for hydrophobic, high blocking feed streams.

Virosart® HC retains  $\geq 4 \log_{10}$  of small non-enveloped viruses (e.g. PPV, MVM) and  $\geq 6 \log_{10}$  of large enveloped viruses (e.g. MuLV). This filter offers excellent virus retention over the entire flow decay profile independently of operating pressure.

Customized process-specific virus LRV determination can be individually provided by our Confidence® Virus Clearance Services (see page 7 for details).

## Integrity Testing

Virosart® HC filters are tested for integrity using a water based diffusion test, e.g. based on the Sartocheck® technology of Sartorius Stedim Biotech. Virosart® HC filters have been validated for logarithmic reduction values of  $\geq 4 \log_{10}$  for small non-enveloped viruses using bacteriophage PP7 as the model virus.

\* Virosart® Max is a specifically optimized virus pre-filter significantly increasing downstream virus filter performance. Virosart® Max is a patented technology (DE 10 2011 105 525 B4) binding aggregates efficiently through hydrophobic interactions with polyamide, independently of process conditions such as conductivity from biological feed streams (mAbs, plasma derivatives or recombinant proteins).

# Technical Data



|  | Minisart®  | Capsule & Midicaps®   |  | T-Style Maxicaps® & Cartridge   |
|--|--|---|--|---|
| <b>Nominal filtration area</b>                           | 5.0 cm <sup>2</sup>  | 240 cm <sup>2</sup><br>0.27 ft <sup>2</sup>   | 0.27 m <sup>2</sup><br>2.9 ft <sup>2</sup> | 0.9 m <sup>2</sup> , 1.8 m <sup>2</sup> , 2.7 m <sup>2</sup><br>9.7 ft <sup>2</sup> , 19.4 ft <sup>2</sup> , 29 ft <sup>2</sup>   |
| <b>To be used for</b>                                    | <ul style="list-style-type: none"> <li>Scale-down work</li> <li>Flow &amp; capacity studies</li> <li>Optimization of pre-filter-   final-filter-ratio</li> <li>GLP spiking studies (IT tested Minisart®)</li> </ul>                      | <ul style="list-style-type: none"> <li>Scale-up studies</li> <li>Small scale production</li> </ul>  |  | <ul style="list-style-type: none"> <li>Large scale manufacturing</li> </ul>   |
| <b>Typical filtration volume</b>                         | < 200 mL   | < 5 L   | < 80 L                                     | > 80 L  |
| <b>Available connectors</b>                              | <ul style="list-style-type: none"> <li>Female luer lock inlet &amp; male luer lock outlet</li> </ul>   | <ul style="list-style-type: none"> <li>Midicaps® &amp; Capsules: FF: 3/4" triclamp (sanitary) connector inlet &amp; outlet</li> </ul>                                     |  | <ul style="list-style-type: none"> <li>Maxicaps®: Sanitary inlet &amp; outlet adapter</li> <li>Cartridges: S-adapter top, 2 flange bayonet adapter with double o-ring bottom</li> </ul>   |
| <b>Sterilization</b>                                     | <ul style="list-style-type: none"> <li>Autoclaving: 121 °C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul> <p>⚠ No inline steaming of Minisart®</p>  | <ul style="list-style-type: none"> <li>Autoclaving: 121 °C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul> <p>⚠ No inline steaming of Capsule &amp; Midicaps®</p> |  | <ul style="list-style-type: none"> <li>Maxicaps®   Cartridges: Autoclaving 121 °C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul> <p>⚠ No inline steaming of Maxicaps®</p> <ul style="list-style-type: none"> <li>Cartridges: Steaming 130 °C @ 2.7 bar   39.6 psi for 60 min up to 3 cycles</li> </ul> |
| <b>Operating parameters</b>                              | <ul style="list-style-type: none"> <li>In the direction of filtration: max. 5.0 bar   73 psi at 20 °C, max. 0.2 bar   2.92 psi at 121 °C</li> <li>In the reversed direction of filtration: max. 0.05 bar   0.725 psi at 20 °C</li> </ul> |   |  |   |
| <b>Water based diffusion test at 4.5 bar   65.25 psi</b> | N/A  | <ul style="list-style-type: none"> <li>2 mL/min (240 cm<sup>2</sup>)</li> <li>8 mL/min (0.27 m<sup>2</sup>)</li> </ul>  |  | <ul style="list-style-type: none"> <li>24 mL/min (0.9 m<sup>2</sup>)</li> <li>48 mL/min (1.8 m<sup>2</sup>)</li> <li>72 mL/min (2.7 m<sup>2</sup>)</li> </ul>   |

# Materials

## Device

### Cartridges, T-Style Maxicaps®, Capsules & Midicaps


|                    |               |
|--------------------|---------------|
| Supportive fleece  | Polypropylene |
| Capsule housing    | Polypropylene |
| End caps           | Polypropylene |
| Core (not capsule) | Polypropylene |

### Membrane

|           |                                   |
|-----------|-----------------------------------|
| Material  | Surface modified polyethersulfone |
| Pore size | 20 nm nominal                     |
| Format    | Double layer                      |

### Minisart®

|                 |               |
|-----------------|---------------|
| Capsule housing | Polycarbonate |
|-----------------|---------------|




Construction of Virosart® HC cartridge and capsule with zoom on cross section and membrane.


# Performance

## Characteristic Water Flow Rates


Virosart® HC 30" Cartridge & 30" T-Style Maxicaps®  
(2.7 m<sup>2</sup> | 29 ft<sup>2</sup>)




Virosart® HC 20" Cartridge & 20" T-Style Maxicaps®  
(1.8 m<sup>2</sup> | 19.4 ft<sup>2</sup>)




Virosart® HC 10" Cartridge & 10" T-Style Maxicaps®  
(0.9 m<sup>2</sup> | 9.7 ft<sup>2</sup>)




Virosart® HC Midicaps® (0.27 m<sup>2</sup> | 2.9 ft<sup>2</sup>)



Virosart® HC Capsule (240 cm<sup>2</sup> | 0.27 ft<sup>2</sup>)



Virosart® HC Minisart® (5 cm<sup>2</sup>)



## Regulatory Compliance

Each individual filter is tested for integrity (except 539VM-----B)

- Validated for  $\geq 4 \log_{10}$  removal of small non-enveloped viruses using bacteriophage PP7
- Designed, developed and manufactured in accordance with an ISO 9001 certified Quality Management System
- Meet or exceed the requirements for WFI quality standards set by the current USP
- Non pyrogenic according to USP Bacterial Endotoxins
- USP Plastic Class Test VI

## Technical References

|                         |   |
|-------------------------|---|
| Validation Guide        | 2838409-000   |
| Extractables Guide      | 2889522-000   |
| Virus Information Guide | SPK5752-e   |
| Application Note        | Autoclaving Virosart® Minisart® devices (SPK4110-e) |

## Ordering Information



### Cartridge

539 25 28 ■

#### Adapter

25: S-adapter top, 2 flange bayonet adapter with double o-ring bottom

#### Filter size

V1: 10" 0.9 m<sup>2</sup> | 9.7 ft<sup>2</sup>  
V2: 20" 1.8 m<sup>2</sup> | 19.4 ft<sup>2</sup>  
V3: 30" 2.7 m<sup>2</sup> | 29 ft<sup>2</sup>



### T-Style Maxicaps®

539 83 28 ■ -- SS

#### Filter size

V1: 10" 0.9 m<sup>2</sup> | 9.7 ft<sup>2</sup>  
V2: 20" 1.8 m<sup>2</sup> | 19.4 ft<sup>2</sup>  
V3: 30" 2.7 m<sup>2</sup> | 29 ft<sup>2</sup>

#### Adapter

SS: Sanitary inlet – and outlet adapter



### Midicaps®

539 53 28 V9 -- FF -- V

#### Filter size

V9: Size 9  
0.27 m<sup>2</sup> | 2.9 ft<sup>2</sup>

#### Adapter

FF: 3/4" triclamp (sanitary) connector inlet & outlet

#### Units per package

V: Two pieces



### Capsule

539 13 28 V4 -- FF -- B

#### Filter size

V4: Size 4  
240 cm<sup>2</sup> | 0.26 ft<sup>2</sup>

#### Adapter

FF: 3/4" triclamp (sanitary) connector inlet & outlet

#### Units per package

B: Five pieces



### Minisart®

539 VM -- -- -- -- ■ ■

#### IT

IT: Integrity tested  
--: Not integrity tested

#### Units per package

A: Four pieces  
B: Five pieces

# Accessories & Services

## Adaptive Pre-Filtration

Virosart® Max\* protects your virus filter irrespective of the process conditions. Virosart® Max will downsize your process and reduce your total virus filtration costs.



## Integrity Testing using Sartocheck®

Fully automated Virosart® integrity testing to guarantee intactness of the Virosart® filter applying pre- and post-use diffusion tests.



## Filter Holders and Housing


Filter holders are designed to accommodate all different Maxicaps® heights. Different standard designs of filter housings are available for cartridges from 10" to 30".

## Single-use Systems

Flexible processing with Flexact® VR system for production from pilot plants up to commercial processing.

## Customized Systems

High level of automation and individual requirements can be realized by customized single-use or hybrid solutions.



## Sartorius Confidence® Virus Clearance Services

Our virus clearance service is the perfect complement to Virosart® HC filters, providing:

- GLP virus clearance studies
- Process feasibility studies
- Regulatory guidance

We use a variety of different relevant and model viruses including BVDV, PRV, HAV and PPV. The combination of product and services provides you with a comprehensive virus clearance solution that gives you the confidence you need to proceed.

## Sartorius Biologics Testing Services

Your partner to assure an effective virus clearance strategy for your blood | plasma product and manufacturing process.

**Germany**

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178

 For more information, visit  
[sartorius.com](https://www.sartorius.com)

Specifications subject to change without notice.

© 2025 Copyright Sartorius Stedim Biotech GmbH, August-Spindler-Strasse 11, 37079 Goettingen, Germany

DIR: 2650742-000-05 | Material-No.: 85037-561-68

Status: 03 | 2025